

STOCK CODE 603387

One Step Test for SARS-CoV-2 Antigen (Colloidal Gold)



CE Marked

Rapid, Reliable, Early and Easy
Detection of COVID-19



Intended Use

The test kit is intended for the qualitative detection of nucleocapsid protein antigen from SARS-CoV-2 in human nasal swab samples or saliva samples from patients suspected of COVID-19 infection by a healthcare provider.




- **PCR-based tests** can detect small amounts of viral genetic materials.
- **Antigen tests** detect the presence of viral proteins and can return positive results when a person is most infectious.
- - - **Antibody tests** detect the body's immune response to the virus.

Characteristics


- Specimen type: nasal swab, saliva
- Rapid test. The test result is available in 10-15 min.
- Read test results visually. Do not require test equipment.
- Early detection of SARS-CoV-2 infection
- Simple operation and highly-efficient test.

Operation

Sample Pretreatment >>







Test >>



⌚
Read result in 10-15 min

Step 2

Test Results




Getein SARS-CoV-2 Testing Products

Nucleic Acid Test



Antigen Test



Quantitative Detection

One Step Test for Novel Coronavirus (2019-nCoV) IgM/IgG antibody (Colloidal Gold)



Qualitative Detection

Antibody Test



Getein Biotech, Inc.

Add: No.9 Bofu Road, Luhe District, Nanjing, 211505, China
 Tel: +86-25-68568508/68568594
 Fax: +86-25-68568500
 E-mail: sales@getein.com.cn; overseas@getein.com.cn
 Web: en.bio-gp.com.cn



Lotus NL B.V.

Add: Koningin Julianaplein 10, 1e Verd, 2595AA, The Hague, Netherlands
 Tel: +31645171879(English), +31626669008(Dutch)
 Email: peter@lotusnl.com

